

Start Up

1. **Enable** the CEE APOGEE Bake Plate on the TNFC Access System.
2. **Turn ON** the machine using the button on the **right**.
3. **Login** using the following credentials Username: **admin** Password: **admin2**
4. **Confirm** that the system is ready for operation by navigating to **Tools>Manual Control**:
 - Plate temperature: ~20 °C • Ambient Temperature: ~22 °C
 - Lift pin height: 5.0 mm • Humidity: ~31-32%
5. **Confirm** that the system is ready for operation by **checking the gas valves**:
 - Main nitrogen gas valve: **Optional (ONLY open when proximity bake is needed)**
 - Hood nitrogen valve: **Optional (ONLY open when hood Nitrogen purge is needed)**
 - Vacuum pump: **Optional (ONLY turn on the pump next to it when vacuum bake is needed)**

WARNING: Bake Plate surface might be **VERY HOT, BE CAUTIOUS**.

REMARKS: The plate temperature at ~20°C is just a reference temperature for the idling state. The process can still start at any plate temperature just that there might be risk of incurring process variation.

REMARKS: The lift pin height (home) can be changed according to user preference. For example, users can set the home height to 0mm through **Manual Control** and the sample will be placed on the bake plate surface directly at the start of the process.

Bake Methods

6. **Vacuum method (Hard contact)**
 - Sample is held securely to the bake plate surface through applying a **vacuum**.
 - **Uniform heating** across the sample surface.
 - Applicable to samples where **back side contact** is **NOT** a problem.
7. **Contact method (Soft contact)**
 - Sample is held to the bake plate surface through **gravity ONLY**.
 - **Less uniform** and **less efficient** heating across the sample surface.
8. **Proximity method (Proximity contact)**
 - Sample **floats** (1-4mm) on the bake plate surface by a **thin nitrogen gas layer**.
 - **Slower warm-up** than contact bake methods

- Advantages when baking thick films where **blistering** would be a problem.
- A high degree of **uniformity** even for **cambered**, or **warped** samples.

Making A New Recipe

9. Navigate to the **Recipes** Tab.
10. Press **New**.
11. **Name** your recipe accordingly.
12. **Set** the **plate temperature** by typing in the input field.
13. **Change** the **time** and **process method** (lift pin height, bake method) accordingly.
14. Insert more steps using the “**Insert**” button if it is necessary.
15. **Save** the recipe after editing.

***REMARKS:** The default recipe editing mode is a relatively “simple” mode that only allows a single temperature target, users need to change to advance editing mode for purposes such as setting up a temperature ramp and holding at a certain temperature for a prolonged period.*

Advance Editing Mode

16. **Click** the “**Advanced**” button, a “!” is going to show up within the button, click again.
17. **Insert** a new step using the “**Insert**” button.
18. **Click** on the newly added step.
19. **Change** the control and the action accordingly.

***REMARKS:** The advance editing mode allows a high degree of freedom on controlling the process. Users can add multiple temperature target within a single process with ramping rate up to **6°C /min**. Users can also add fixed time delays to make the machine hold at a certain temperature according to their need. Changing pin height during the process is also possible.*

20. **Press** “**Update**” to log the change and a new step will be created.
21. **Organize** the order of the steps using the “**up**” and “**down**” arrows.
22. **Modify** the tolerances and preconditions using their respective button if necessary.
23. **Save** the recipe after editing.

Loading And Running Existing Recipe

24. **Navigate** to **Recipes** Tab.

25. Press “Load”.
26. Click on the designated **recipe name**.
27. Press “Run”.
28. Start the process using the “**Start**” button on the main **Process Tab**.

Loading A Sample

29. **Open** the hood using the handle.
30. **Align** and **place** the sample on the lift pin.
31. **Close** the hood **slowly** and **gently**.
32. **Turn on** the hood nitrogen valve if **nitrogen purging** is necessary.

Sample Unloading

33. Once the process is **completed**, the machine will make a **beeping noise**.
34. *Check the **bake plate temperature** from **Tools>Manual Control**.*
35. *Turn off the hood nitrogen purging valve if it was opened.*
36. *Open the hood **ONLY** when it is cooled down to a **safe temperature**.*

WARNING: The **hood handle** might be **HOT** when the interior temperature is **HIGH**, **BE CAUTIOUS**.

REMARKS: It may take some time for it to cool down, please **BE PATIENT**.

37. *Remove the sample from the bake plate.*

Shut Down

38. **Confirm** that the gas valves & vacuum are off:
 - Main nitrogen gas valve: **OFF**
 - Hood nitrogen valve: **OFF**
 - Vacuum pump: **OFF**
39. **Logout** from the bake plate by **admin>Log Out**.
40. **Turn OFF** the machine using the **right** button.
41. **Disable** the **CEE Apogee Bake Plate** on the **TNFC Access System**.